

## PH #1009 Point-of-Use Vial Filling

### Application: VIAL FILLING APPLICATION

Many Pharmaceutical manufacturers have requirements for fine temperature control at point-of-use. This customer required precise temperature control of the product during vial filling. They require double tubesheet models and full documentation packages.

### Solution

- Sanitary Shell & Tube Heat Exchanger
- Double Tubesheet, 73 Series
- ASME BPE sanitary flanges



### Features & Benefits

- DTS eliminates risk of cross contamination
- No dead legs or crevices
- Fully drainable
- Delivered from stock



#### How to contact Exergy

Please call 1-516-832-9300  
or visit [www.exergyllc.com](http://www.exergyllc.com)  
to locate your local rep.